

Data Injection Checker

User Guide - EXE Version

Version 1.0

Purpose

DataInjectionChecker.exe validates that data from RAW interferent result files has been correctly injected into the reference workbook. The application automatically compares RAW files with the Data Entry worksheet and identifies missing, extra or inconsistent values.

No Python Installation Required

End users only need to launch DataInjectionChecker.exe. Python and external libraries are not required on the target workstation.

Folder Preparation

Place DataInjectionChecker.exe, one reference workbook named DataAnalysisSheet_*.xlsx and all RAW Excel files in the same folder.

Reference Workbook Detection

The tool automatically detects a reference file matching DataAnalysisSheet_*.xlsx. Only one reference workbook should be present in the folder.

Supported Interferents

Protein, Triglycerides (TG), Hemoglobin (HB), Rheumatoid Factor (RF), Bilirubin unconjugated (BR) and Bilirubin conjugated (BR CONJ).

Running the Tool

Double-click DataInjectionChecker.exe. The application automatically reads the RAW files, compares them against the reference workbook and generates a validation report.

Generated Report

The application creates Injection_Check_Report.xlsx. If a report already exists, a new file is created automatically using an incremental numbering scheme (Injection_Check_Report_01.xlsx, Injection_Check_Report_02.xlsx, etc.). Existing reports are never overwritten.

Report Worksheets

Summary: overall validation result. Detailed Comparison: file statistics. Missing Triplets: values expected but not found. Extra Triplets: unexpected values found in the reference workbook.

Status Definitions

■ OK: all values match. ■ CHECK REQUIRED: differences were detected. ■ ERROR: the file could not be processed.

Recommended Workflow

Launch the application, open the generated report, review the Summary worksheet and investigate any CHECK REQUIRED result using the Missing Triplets and Extra Triplets worksheets.